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Léo Duvivier  
Inside Sales Specialist - Europe

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
Throughout Users, Labs, and Sites

## Traceability & Compliance

Costly Errors and Usability Challenges in GxP/Regulated Environment

## Productivity & Efficiency

Effectiveness of the Laboratory

## Ergonomics

Incidences of Repetitive Strain Injury (RSI)

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
- 0.2-10 µL
- 5-120 µL
- 10-300 µL
- 50-1,200 µL

## Modes de pipetage

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Pipetting	✓
Reverse Pipetting	✓
Manual Pipetting	✓
Multi Dispensing	✓
Diluting	✓
Sequential Dispensing	✓
Multi Aspiration	✓
Titrate	✓

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**OneLab**   
design & execute

Build A  
System  
For Your  
Application

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## Peltier+

### Fast Heating and Cooling

Rapid cooling down to 0°C and heating up to 99°C

## Vacuum+

### Automated Solid Phase Extraction

Enhanced vacuum control designed for use with Waters uElution plates

## Magnet+

### Simple, Reliable Separation

Avoid centrifugation, a bottleneck in many automated processes

## Shaker+

### Rapid, Precise Mixing

Vibration free motion ensures highest mixing speeds on market  
200-3000 rpm

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the pipetting robot

## Examples

AccQ-Tag Ultra LC **Amino Acid** Analysis

Immunosuppressant sample prep


Associates of Cape Cod, Inc.

Fungitell® **Glucan** Detection Assay

Sample Preparation of mAbs  
for **Middle-level** Analysis

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

# OneLab Evolving Connected Environment

**Andrew**  
Liquid Handling Robot

**Ergonomic Bluetooth  
Electronic Pipettes**

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


ALL TUBES & VIALS MICROPLATES BOTTLES & RESERVOIRS MICROCHIPS & DEVICES COLUMNS

Filters

- Number of wells
  - 1 6 8
  - 12 24 36
  - 48 54 60
  - 96 384
- Well capacity
  - < 100 uL 100 uL - 1 mL
  - 1 mL - 50 mL
  - 50 mL - 250 mL > 250 mL
- Device compatibility
  - Absorbance plate reader
  - Andrew+
  - Heater-Shaker+
  - Magnet+
  - Peltier+
  - Shaker+
  - Vacuum+
  - PCR thermal cycler
  - Photoredox

Frequently used items

- Fisherbrand™ Premium 1.5 mL microtube [Add](#)
- Falcon® 50 mL conical centrifuge tube [Add](#)
- Falcon® 15 mL conical centrifuge tube [Add](#)
- Fisherbrand™ 2 mL skirted conical microtube [Add](#)
- 10 mL flat-bottom screw cap test tube [Add](#)
- 13x84mm collection tube, virus preservation... [Add](#)
- 16x100mm collection tube, COPAN UTM®... [Add](#)
- 16x100mm collection tube, IMPROVIRAL™... [Add](#)
- 16x100mm collection tube, virus preservation... [Add](#)
- 2 mL flint glass vial [Add](#)
- ACC, Pyrotube® 16x90mm depyrogenated... [Add](#)
- Agilent 12x32mm screw-top vial, 0.3 mL insert in Water... [Add](#)
- Agilent 2 mL 12x32mm crimp-top vial [Add](#)
- Agilent 2 mL 12x32mm crimp-top vial in Genevac 48... [Add](#)
- Agilent 2 mL 12x32mm crimp-top vial in Waters 48x... [Add](#)
- Agilent 2 mL 12x32mm screw-top vial in 54x vial plate [Add](#)
- Agilent 2 mL 12x32mm screw-top vial in Vanquish™... [Add](#)
- Axygen® 1.5 mL Maxymum Recovery® snaplo... [Add](#)

3 labware to add on the bench

- Fisherbrand™ Premium 1.5 mL microtube [Add](#)
- Eppendorf twin.tec® 96-well skirted LoBind®... [Add](#)
- INTEGRA 10 mL multichannel reservoir [Add](#)

Fisherbrand™ Premium 1.5 mL microtube

Labware representation

[Add to labware list](#) [Manufacturer data](#)

Labware part number	11926955
Max volume	1.6 mL
Working volume	1.5 mL
Dead volume	20 µL
Bottom shape	V-shaped
Labware manufacturer	Fisher Scientific
Closure	Snap cap
Sterility requirements	Autoclavable, DNase-Free, RNase-Free, Pyrogen-Free
Material	PP
Surface treatment	Highly polished
Color	Clear

Compatible labware holders

- Microtube 218.2152 [Available in this lab](#)
- 1.5mL Microtube Cooled 218.2451

Compatible pipettors

1-channel pipettes

Channel / volume range (µL)	Requirements
1-ch 0.2-10 µL SKU: 186009769	Fully compatible <a href="#">Available in this lab</a> <a href="#">Shop on waters.com</a>

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Kairos Derivatization, 48 Samples - Andrew+ v. 1 Execute

AccQ-Tag 3x Label

Eppendorf twin.tec® 96-well plate

Waters 2 mL 96-well collection plate

Agilent 6-column reagent reservoir

INTEGRA 10 mL multichannel reservoir

STEPS LABWARE INFO

1700 rpm after 10 seconds

20. Dispense 1 mL from AccQ-Tag 3x Label to INTEGRA 10 mL multichannel reservoir

21. Dispense 500 µL from AccQ-Tag 3x Label to INTEGRA 10 mL multichannel reservoir

Leftover of 0 µL in source AccQ-Tag 3x Label (deadvolume: 20 µL).

22. Dispense 20 µL from INTEGRA 10 mL multichannel reservoir to Eppendorf twin.tec® 96-well plate A1:H1

23. Dispense 20 µL from INTEGRA 10 mL multichannel reservoir to Eppendorf twin.tec® 96-well plate A2:H2

24. Dispense 20 µL from INTEGRA 10 mL multichannel reservoir to Eppendorf twin.tec® 96-well plate A3:H3

25. Dispense 20 µL from

# La creation des protocoles via OneLab

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### Pipetting mode

Forward  Reverse  Repetitive mode ⓘ

Blow-out ⓘ

Pause before blow-out

### Handling liquid viscosity

Pipetting speed fine adjustment

### Aspiration speed

Slow  Normal  Fast

### Dispensing speed

Slow  Normal  Fast

### Air cushion

No air cushion

Air top cushion (high viscosity aspiration) ⓘ

Air bottom cushion (low viscosity aspiration) ⓘ

### Pipette moving speed ⓘ

Slow  Normal

### Tip position ⓘ

#### Sources

With respect to liquid  
 With respect to bottom  
 Avoid touching bottom

#### Destinations

On-the-fly  
 With respect to liquid  
 With respect to bottom  
 Avoid touching bottom

### Tip choice ⓘ

Change tip at the beginning

Change tip between pipetting steps

Use tips with filter

### Mixing ⓘ

Sources

#### Mixing

[-] 3 times [+]

#### Speed

Slow  Normal  Fast

#### Volume

0  $\mu$ L

Destinations

#### Mixing

[-] 3 times [+]

#### Speed

Slow  Normal  Fast

#### Volume

0  $\mu$ L

### Automation

#### Mechanical arm motion speed

Normal  Fast

#### Custom tip position ⓘ

##### Sources *With respect to liquid*

Customized aspiration height

[-] 0 mm [+]

##### Destinations *With respect to liquid*

Customized dispensing height

[-] 0 mm [+]

#### Tip geometry control

Verify each tip used ⓘ

### Guidelines

Normal : **B** *I* U ☰ ☷ ☹ ☎ 📷 📅

Provide additional guidance about the execution of this step, e.e. how to proceed etc.

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
### Actions

Select an action to apply on labware

 Pipetting	 Serial dilution	 Concentration normalization
 Centrifugation	 Apply magnetic field	 Flow-through extraction
 Vortexing	 Shaking	 Heat / cool
 Discard Solution	 Timer	 User action / notification

### 2. Serial dilution

Drag required labware for serial dilution Continue

Patient 1 Patient 2 Patient 3

Patient 4 Patient 5 Patient 6

Patient 7 Patient 8

	1	2	3	4	5	6	7	8	9	10	11	12
A	●	○	○	○	○	○	○	○	○	○	○	○
B	●	○	○	○	○	○	○	○	○	○	○	○
C	●	○	○	○	○	○	○	○	○	○	○	○
D	●	○	○	○	○	○	○	○	○	○	○	○
E	●	○	○	○	○	○	○	○	○	○	○	○
F	●	○	○	○	○	○	○	○	○	○	○	○
G	●	○	○	○	○	○	○	○	○	○	○	○
H	●	○	○	○	○	○	○	○	○	○	○	○


Buffer

You want to dilute: ⓘ

One sample

A column of 8 samples

— Samples to dilute

Select wells

Drop one compatible labware

— Dilution buffer

Buffer  
AUTO (250 µL)

Drop one compatible labware

— Destinations

Select wells

Drop one compatible labware

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 Edit protocol  Step list  Cancel

### Select the number of samples

Minimum 1 sample, maximum 8 samples

The maximum number is based on the sample references listed in the protocol

Adapt for 8 samples

or

### Import a protocol worklist file (.csv)

 **Help:** Prepare your protocol worklist for execution

[Download template](#) to create a compatible import file

Drag and drop a csv file

Select a csv file

Alternatively, you can execute the protocol with the existing parameters

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

**Lab productivity**

**Reproducibility** between assays

**Comparability** between manual and automated preparations

**Reliable** sample preparation

**Simple & easy** automation

Minimizes **error!**

Walk-away **time!**

Waters™