

## PT 86b: Waters for hospital use – Microbiology of dialysis waters (hemodiafiltration)

### Details of the program:

- Proficiency testing scheme created in 2019
- 2 rounds per annual series
- The time for analysis is 3 weeks



### Schedule:

CODE	MATRIX	WEIGHT	NEXT SHIPMENT
0286	Dialysis water* ( <i>ultra-pure water similar to dialysate</i> )	120 ml	2024/01/15
0286	Dialysis water* ( <i>ultra-pure water similar to dialysate</i> )	120 ml	2024/05/27

\* Refrigerated samples

### ANALYTES

Total viable count at 22°C for 7 days on medium similar to R2A media  
Identification of indicator microorganisms

### METHODS

Membrane filtration and biochemical or genomic test (identification)

*Note: Matrices and analytes may be changed or removed for technical or scientific reasons.  
Please refer to current application form available in your member area ([www.bipea.org](http://www.bipea.org)).*