

PT 105a: Microbiology in feed – *Salmonella* spp detection

Details of the program:

- Proficiency testing scheme created in 2020
- **15 registered laboratories from 10 countries**
- 2 rounds per annual series
- The time for analysis is 3 weeks
- Samples are shipped via express carrier.



Schedule:

CODE	MATRIX	VOLUME	NEXT SHIPMENT
01105	Feed – <i>Salmonella</i> *	25g	2023/10/23
02105	Feed – <i>Salmonella</i> *	25g	
03105	Feed – <i>Salmonella</i> *	25g	
01105	Feed – <i>Salmonella</i> *	25g	2024/05/06
02105	Feed – <i>Salmonella</i> *	25g	
03105	Feed – <i>Salmonella</i> *	25g	

*Refrigerated samples

ANALYTES

Salmonella spp detection

Note: Matrices and analytes may be changed or removed for technical or scientific reasons.
Please refer to current application form available in your member area (www.bipea.org).