

PT 71: Allergens in wines

Details of the program:

- Proficiency testing scheme created in 2014
- **31 registered laboratories from 10 countries**
- **PTS accredited by COFRAC**
- 3 rounds per annual series
- The time for analysis is 3 weeks
- Samples are shipped via express carrier.



Schedule:

SHIPMENT DATE	CODE	MATRIX	VOLUME
2022/10/03	1171	Rosé wine spiked with Ovalbumine*	250 ml
	1271	Rosé wine spiked with Casein*	250 ml
	1371	Rosé wine spiked with Lysozyme*	250 ml
	1471	Rosé wine spiked with Beta-lactoglobuline* NA	250 ml
	1771	Rosé wine spiked with Gluten* ^{NA}	250 ml
2023/02/13	0571	White wine spiked with Ovalbumine*	250 ml
	0671	White wine spiked with Casein*	250 ml
	0771	White wine spiked with Lysozyme*	250 ml
	1571	White wine spiked with Beta-lactoglobuline* NA	250 ml
	1871	White wine spiked with Gluten* ^{NA}	250 ml
2023/05/22	0871	Red wine spiked with Ovalbumine*	250 ml
	0971	Red wine spiked with Casein*	250 ml
	1071	Red wine spiked with Lysozyme*	250 ml
	1671	Red wine spiked with Beta-lactoglobuline* NA	250 ml
	1971	Red wine spiked with Gluten* ^{NA}	250 ml

* Refrigerated samples

ANALYTES

Casein
 Lysozyme
 Ovalbumine
 Egg proteins^{NA}
 Beta-lactoglobuline^{NA}
 Gluten^{NA}

NA *: Not accredited parameter

Note: Matrices and analytes may be changed or removed for technical or scientific reasons.

Please refer to current application form available in your member area (www.bipea.org).