

## PT 74: Melamine and cyanuric acid in Milk

### Details of the program:

- Proficiency testing scheme created in 2015
- **24 registered laboratories from 12 countries**
- **PTS accredited by COFRAC**
- 2 rounds per annual series
- The time for analysis is 4 weeks
- Samples are shipped via express carrier at the end of the previous month.



### Schedule:

DATE	CODE	MATRIX	WEIGHT
September	0174	Milk powder	150 g
March	0174	Milk powder	150 g

### ANALYTES

Range: 0.1 to 2.5 mg/kg

Melamine  
Cyanuric acid

*Note: Matrices and analytes may be changed or removed for technical or scientific reasons.  
Please refer to current application form available in your member area ([www.bipea.org](http://www.bipea.org)).*

